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The screenshot displays the Mosaic web interface for 'Order 192' (Kinase 2 RJF) by user 'Fry, Richard'. The interface is divided into a left-hand navigation menu and a main content area. The navigation menu includes sections for Orders, Fulfilment, Inventory, Lists/Sets, Reporting, Admin, and Help. The main content area features tabs for Status, Streams, Compounds, and Options. The 'Streams' tab is active, showing a table of assay streams with columns for '#', 'Summary', and 'Can have child streams'. Three streams are listed, each with a checkbox and a plus icon. Below the table are 'Move Up' and 'Move Down' buttons. A detailed configuration panel for 'Kinase IC50 10pt 10mM' is shown, with fields for Parent (0), Number of Copies (1), Destination Type (Despatch), Recipient (Fry, Richard), Delivery Location, Target Amount (20 µL), Top Concentration (10 mM), Solvent (100% DMSO), Plate Type (96 well, shallow, v-bottom polystyrene plate), Layout (10pt IC50 NC - 8 Compounds per Plate), Dilution Series (1:2 (10 doses)), Target (Kinase Panel), and Project (Kinase Panel IC50). An 'Apply Changes' button is at the bottom.

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1. ["High Throughput Screening , Methods and Protocols " William P Janzen. Humana Press ISBN 0-89603-889-0](#)
2. Overview of Liquid Handling Instrumentation for High-Throughput Screening Applications, Stewart Rudnicki, Sean Johnston, *Curr. Protoc. Chem Biol.* 1:43-54. © 2009 by John Wiley & Sons, Inc
3. Lipinski, C (2005). Compound precipitation from DMSO and the synergy between water uptake and freeze/thaw cycles. Oral Presentation at LRIG Meeting, January 2005
4. Harrison, WJ "The importance of automated sample management systems in realising the potential of large compound libraries in drug discovery". *J. Biomol. Screen* 2.203
5. Report: "Small-Scale Benchtop Automation Trends 2011" Dr. John W. Comley, HTSTech Ltd

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